

WHEAT PACKS

HEAVEN SCENT

Therapeutic Hot & Cold Packs

- Naturally Treated Export quality Grain
- Available with or without Soothing Lavender Scent
- Natural Hot & Cold Relief for:
 - Headaches, Aching joints, Back & Arthritic pain & more



New!

Online ordering!

You can prepare
and email an order
from this catalogue

Details next page...

ARCHERS
Healthcare

Fax 1300 78 66 00
www.archers.net.au

To order online...

1. Type your details below
2. Put quantities in the QTY boxes next to the products you want.
3. Click the Send Email button on any page when you've finished.
(Only do this once, it attaches this form to an email addressed to us)
4. Send your email.

Pharmacy Name:

Suburb:

Wholesaler:

Account:

Contact:

* All fields highlighted
In red are required.

Heaven Scent Hot & Cold Wheat Packs



HSLP Large Wheat Pk Lavender

HSLPLF Large Wheat Pk Lavender Free

(21 cm x 16cm)

HSSP Small Wheat Pk Lavender

HSSPLF Small Wheat Pk Lavender Free

(38cm x 16cm)

HSNW Neck Wheat Pk Lavender

HSNWLF Neck Wheat Pk Lavender Free

(50cm x 16cm)

HSEP Eye Wheat Pk Lavender

HSEPLF Eye Wheat Pk Lavender Free

(19cm x 10cm)

Illustrations are examples only.

Actual print patterns may vary with supply.

IMPORTANT HEATING INSTRUCTIONS

The following is a recommended guide for safely heating your Heat Pack.

1. Your Heat Pack should only be heated in a Microwave oven equipped with a revolving turn table.
2. Never overheat your Heat Pack. It should never be too hot and uncomfortable to handle.
3. Because Microwave ovens vary in power output, we strongly suggest you first start with **TWO minutes**. Check your **Heat Pack** after this time and if further heating is required increase by only 15 Seconds at a time until desired temperature is achieved. Do not overheat. If you are aware of a burnt odour, this indicates overheating and the pack should be replaced.
4. We suggest once you have found the correct healing time, your **Heat Pack** can also be used to warm those cold feet, cold hands and back.
5. Do not wrap your heat Pack in towelling or the like when applying to the body. For best results, just apply it as it is.
6. Do not re-heat the pack until completely cool.
7. Your **Heat Pack** is not suitable for young children (generally under 4 years old), the infirmed, unconscious patients or people who are unable to react to extremes of temperature.
8. With constant use the natural moisture in the grain may dry out. We suggest occasionally placing a glass of water with the pack when heating. This will prolong the life of your Heat Pack.
9. **CLEANING**. Spot clean only, do not immerse in water.

